# Clinpro

Sellante de resina de fotocurado

## Descripción

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### Etch Enamel:

- Using syringe tip, or fiber tip, apply a generous amount of etchant to all enamel surfaces to be sealed, extending beyond the anticipated margin of the sealant.
- Etch for a minimum of 15 seconds, but no longer than 60 seconds.


### Rinse Etched Enamel:

- Thoroughly rinse teeth with air/water spray to remove etchant.
- Do not allow patient to swallow or rinse. If saliva contacts the etched surfaces, re-etch for 5 seconds and rinse.

### Dry Etched Enamel:

- Thoroughly dry the etched surfaces.
- Air should be oil and water free.
- The dry etched surfaces should appear as a matte frosty white. If not, repeat steps 1 and 2. Do not allow the etched surface to be contaminated.

### Apply Sealant:

- Using the syringe needle tip, or a brush, apply sealant into the pits and fissures. Do not let sealant flow beyond the etched surfaces.
- Stirring the sealant with the syringe-tip during or after placement will help eliminate any possible bubbles, and enhance the flow into the pit and fissures. An explorer may also be used.

### Light-Cure:

- Cure the sealant by exposing it to light from a 3M™ ESPE™ Light-Curing Unit, or other curing unit of comparable intensity.
- A 20-second exposure is needed for each surface. The tip of the light should be held as close as possible to the sealant, without actually touching the sealant.
- When set, the sealant forms a hard, opaque film, light yellow in color with a slight surface inhibition.

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